



Indigo

PTA Rapid Exchange Balloon Catheter



Indigo RX - PTA Rapid Exchange Balloon Catheter is a single use, double lumen dilation balloon catheter for peripheral interventions. The device incorporates the flexibility of rapid exchange design with high burst balloons, multiple working lengths and balloon sizes to facilitate usage in most peripheral intervention requirements.

Indigo OTW - Over-the-Wire Balloon Catheter incorporates the rigidity of over-the-wire design and high pressure balloon to allow clinicians a device for optimal balloon performance during peripheral interventions.



Technical Specifications		
	Indigo Rx	Indigo OTW
Balloon Material	Nylon	Nylon
Balloon Material Type	High pressure Semi compliant	High pressure Semi compliant
Balloon Folding	Tri fold	Tri / Quad fold
Minimum Guiding Catheter	5Fr. for 4.5 mm, 5 mm Φ 6Fr. for 6.0 mm, 10 mm Φ	5Fr for 4.5 mm, 5 mm Φ 6Fr for 6.0mm, 10mm Φ
Maximum Guidewire	0.018"	0.035"
Radiopaque Markers	2	2
Working lengths	40 cm, 60 cm, 80 cm, 130 cm	40 cm, 60 cm, 80 cm, 130 cm
Size - Length	20, 40, 60, 80 mm	20, 40, 60, 80mm
Diameter (Φ)	4.5, 5.0, 6.0, 7.0 mm	4.5, 5.0, 6.0, 7.0, 8.0, 9.0, 10.0 mm

CE marked
Available in custom sizes, and on OEM and OBL basis